


Description

- Sample Tubes with External Threads 5ml 12x92 mm

Part number

- 040851

Pieces by sales : 1000



Description

- Designed for the storage and transportation of biological material. Manufactured from non-toxic polypropylene, the tube provides strength and clarity and exhibits some unique design features. The vial has external threads, providing a smooth and uniform inner surface, thus reducing the risk of contamination. Tubes can be autoclaved (121 °C) in an upright position with caps loosened.
- Made of polypropylene
- For storing and transporting biological materials
- Withstands temperatures from -196°C to 121°C
- Autoclavable
- Non-sterile
- Round bottom tubes only can be centrifuged up to 17,000 x g.
- With white writing area.
- Autoclavable (121 ° C) en position verticale avec les bouchons desserrés.
- Volume : 5 mm
- Size : 12.5 x 93 mm