

**Product Name:** EasyDip™

**Catalogue No.:** M906-12AS

This document replaces any previous version

## 1. Product Description:

- Slide Staining Jars with staining rack and aluminum holder
- Kit includes:
  - 2 white, 1 Blue, 1 Green, 1 Pink and 1 Yellow Jars (M900-12 series)
  - 1 Staining Rack (M905-12DGY)
  - 1 aluminum rack (M906)

## 2. Packaging:

- Case: 1 kit

## 3. Product Specifications:

- Material:
  - Acetal jars
  - Aluminum rack
- Temperature range: -80°C to +30°C
- Each jar can hold up to 12 slides with the staining rack

## 4. Standards and Conformity:

- ISO 2859-1: Sampling and inspection procedures
- CONEG / RoHS: Plastics and colorants are in conformity with CONEG / RoHS standards for heavy metals
- REACH (SVHC): Plastic is in conformity to REACH standards
- LATEX: Material is Latex Free

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic
- Visual attributes
- Dimension measurements

**6. Traceability:**

- **Lot No. Composition:** 8 or 9 digits
  
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label
  2. On label inserted inside the master case
  3. On the inner bag

**7. Storage Conditions:**

- Store at room temperature in normal warehouse conditions
- Protect from any damage to the packaging

**8. Recommended Use:**

- Avoid picric acid and phenol: prolonged use in these reagents will cause degradation
- Follow manufacturer's instructions when using colorant or other chemical substance

**OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*