



Certified Quality
Management System

GOSSELIN™
A Corning Brand

TECHNICAL DATA SHEET

**DESCRIPTION : STRAW PIPETTE 2.2ML PP NA L240,25 P
ER BAG, STERILE SAL 10-3**

Caractéristiques techniques

- **Nominal capacity**
2,2 ml (l 240 mm)
- **Nominal Capacity (Minimum)**
N/A
- **Nominal Capacity (Maximum)**
N/A
- **Pipette OD**
N/A
- **Tip ID**
N/A
- **Mouthpiece OD**
N/A
- **Plug color**
N/A
- **Graduations**
N/A
- **Length (mm)**
240.00
- **Diameter**
4.50
- **Full Capacity**
2,72 ml
- **Physical resistance**
Flexible to rigid, good mechanical resistance, high chemical resistance
- **Use temperature main body**
-20°C/+120°C

Conditionnement

- **Item(s)/bag**
N/A
- **Storage Condition (°C)**
Room temperature (from 15 to 30°C)
- **Quantity per carton**
6000
- **Packaging**
Bag of 25 items

Références normatives

■ **Non-haemolytic**

N/A

■ **Non-pyrogenic**

N/A

■ **Non-cytotoxic**

N/A

■ **Sterility Assurance Level based on ISO 11137**

N/A

Précaution d'utilisation

- **Single use**
-

Domaines d'application

- **Industries**

- **Environment**

- **Microbiology**

- **Dairy**
-

Disponibilité produit

- **Europe Middle East Africa**

- **North America**

- **Asia Pacific**

- **Latin America**

DOMINIQUE DUTSCHER SAS