

CORNING	Life Sciences	Document Number: LSR00192	Rev.: 1
Description (Class): Reference, Customer Technical Data Sheet – Serological Pipets, 10 mL 0.01 mL Increments			Page: 1 of 2

APPLICATIONS:

For tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in an aseptic environment.

FEATURES:

- Consistent performance
 - Reliability, accuracy is +/-2% at full volume
- DNase/RNase Free:
 - Tested by nuclease assay method and found to be free of detectable DNase/RNase contamination. The assay detection limit is 10^{-7} Kunitz units/uL for DNase and 10^{-9} Kunitz units/uL for RNase.
- Non-pyrogenic
 - Low endotoxins
- Color coding
 - Ease of identification on packaging and blow-out rings
- Reverse graduations
 - Easy reading of the amount of liquid removed or the amount remaining
- Negative graduations
 - Can measure past the stated capacity
- Crisp, legible, rotogravure printing
 - Easily read sample volumes
- Crystal grade polystyrene
 - Can clearly see fluid sample
- Individually wrapped packaging:
 - Two ways to open: Peel or Pop-through

REGULATORY COMPLIANCE:

These products are manufactured under the current ISO 9001 and ISO 13485 Standards, and the current FDA Quality System Regulation 21 CFR Section 820.

STERILITY:

Product is gamma irradiated

CORNING	Life Sciences	Document Number: LSR00192	Rev.: 1
Description (Class): Reference, Customer Technical Data Sheet – Serological Pipets, 10 mL 0.01 mL Increments			Page: 2 of 2

10 mL, 0.01 mL Increments, Serological Pipet

Material: Pipet - Polystyrene; Plug - Polyester fiber

Related Catalog Numbers:

Sterile
Sterile

Cat. No.

357773
357779

Packout Description

50/bag, 200/case
25/bag, 500/case

Dimensions: = XXX mm (XXX inches)

TIP:

Length: = 37.36 (1.471)
Interior Diameter: = 1.57 (.062)
Outside Diameter: = 3.35 (.132)

MOUTHPIECE:

Length: = 25.40 (1.00)
Interior Diameter: = 5.72 (.225)
Outside Diameter: = 8.05 (.317)
Plug Length = 15.88 (.625)

EXTRUSION:

Interior Diameter: = 7.95 (.313)
Outside Diameter: = 9.60 (.378)
Length = 285.75 (11.250)

PIPET SIZE:

10 mL

INCREMENTS:

0.10

GRADUATIONS:

Reverse

OVERFILL CAPACITY:


3.00 mL

COLOR CODE:

Red

GENERAL DIMENSIONS:

A = 346.81(13.65)
B = 3.30 (.130) Diameter
C = 8.05 (.317) Diameter
D = 9.61 (.370) Diameter



Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures.

CORNING	Life Sciences	Document Number: LSR00192	Rev.: 1
Description (Class): Reference, Customer Technical Data Sheet – Serological Pipets, 10 mL 0.01 mL Increments			Page: 1 of 1

Revision History				
Rev	CC#	Description of Change	Name	Date
1	CC-02578	Initial Release	A Burgess	7/29/2014