

CORNING	Life Sciences	Document Number: LSR00204	Rev.: 1
Description (Class): Customer Technical Data Sheet – 2.2 mL Pipet			Page: 1 of 2

APPLICATIONS:

For tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in a sterile environment.

FEATURES:

- Consistent performance
 - Reliability, accuracy is +/-3% at full volume
- DNase/RNase Free:
 - Tested by nuclease assay method and found to be free of detectable DNase/RNase contamination. The assay detection limit is 10^{-7} Kunitz units/uL for DNase and 10^{-9} Kunitz units/uL for RNase.
- Non-pyrogenic
 - Low endotoxins
- Color coding
 - Ease of identification on packaging and blow-out rings
- Crisp and legible printing
 - Easily read sample volumes
- Crystal grade polystyrene
 - Can clearly see fluid sample
- Individually wrapped packaging
 - Two ways to open: Peel or Pop-through

REGULATORY COMPLIANCE:

These products are manufactured under the current ISO 9001 and ISO 13485 Standards, and the current FDA Quality System Regulation 21 CFR Section 820.

STERILITY:

Product is gamma irradiated

Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit www.corning.com/lifesciences/trademarks.

© 2013 Corning Incorporated

All other trademarks in this document are the property of their respective owners.

CORNING	Life Sciences	Document Number: LSR00204	Rev.: 1
Description (Class): Customer Technical Data Sheet – 2.2 mL Pipet			Page: 2 of 2

2.2 mL Pipet

Material: Pipet - Polystyrene; Plug - Polyester fiber

Related Catalog Numbers:

Sterile


Cat. No.

357785

Packout Description

25/bag, 1000/case

DIMENSIONS = XXX mm (XXX inches)			
TIP:		PIPET SIZE:	
Length:	= N/A	2.2 mL	
Interior Diameter:	= 1.63 (.064)	INCREMENTS:	
Outside Diameter:	= 3.73 (.147)	N/A	
MOUTHPIECE:		GRADUATIONS:	
Length:	= N/A	1.0, 2.0, 2.1, 2.2 mL	
Interior Diameter:	= 4.17 (.164)	OVERFILL CAPACITY:	
Outside Diameter:	= 6.27 (.247)	N/A	
Plug Length:	= N/A	COLOR CODE:	
EXTRUSION:		GENERAL DIMENSIONS:	
Interior Diameter:	= N/A	A = 282.58 (11.125)	
Outside Diameter:	= N/A	B = 3.73 (.147) Diameter	
Length:	= N/A	C = 6.27 (.247) Diameter	
		D = 6.27 (.247) Diameter	



Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit www.corning.com/lifesciences/trademarks.

© 2013 Corning Incorporated

All other trademarks in this document are the property of their respective owners.

CORNING	Life Sciences	Document Number: LSR00204	Rev.: 1
Description (Class): Customer Technical Data Sheet – 2.2 mL Pipet			Page: 1 of 1

REVISION HISTORY				
Revision	Date	Description	CC No.	Name
1	11/16/17	Original issuance	CC-07917	C. Roberts

Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.