

TECHNICAL DATA SHEET

Part Number	041513
Description	Thermo Scientific™ ART™ Barrier Pipette Tips in Hinged Racks



Benefits :

- ✓ Perform pipetting accurately and consistently with these barrier pipette tips in hinged racks, which seal when exposed to potential contaminants. ART hinged racks are available in low retention, micropoint, extended length and other specialty styles.
- ✓ Dimensions : 50,7 mm
- ✓ Volume : 100 µL
- ✓ Low retention

Sales unit :

960 Tips
