

## TECHNICAL DATA SHEET

<b>Part Number</b>	041534
<b>Description</b>	Thermo Scientific™ ART™ Barrier Pipette Tips low retention



### Benefits :

- ✓ Unlike standard filter tips, Thermo Scientific™ ART™ Barrier Hinged Rack Pipette Tips seal when exposed to potential contaminants. The pipette tips from the reload inserts fit securely into the environmentally friendly, reusable hinged racks and the covers may be used as hinged or lift-off style. The locking clasp prevents contamination by keeping remaining tips sterile and offers an easy, one-handed operation.
- ✓ Features: Stackable hinged racks save space in your lab
- ✓ Sterile tips are certified free from RNase, DNase and pyrogen
- ✓ Recommended for: Liquid Handling, Molecular Biology, PCR, Biotechnology
- ✓ Sterile
- ✓ Volume : 10 µL
- ✓ Tip 0.5-20 µL,
- ✓ Dimensions : 50.7 mm
- ✓ Low retention

### Sales unit :

960 Tips