

TECHNICAL DATA SHEET

Part Number	041626
Description	ART™ Barrier Extended Length Pipette Tips



Benefits :

- ✓ Unlike standard filter tips, Thermo Scientific™ ART™ Barrier Extended Length Pipette Tips in Hinged Racks seal when exposed to potential contaminants. The pipette tips from the reload inserts fit securely into the environmentally friendly, reusable hinged racks and the covers may be used as hinged or lift-off style. The locking clasp prevents contamination by keeping remaining tips sterile and offers an easy, one-handed operation. ART hinged racks are available in low retention, micropoint, extended length, and other specialty styles.
- ✓ Sterile
- ✓ Dimensions : 11.85 cm
- ✓ Volume : 1000 µL
- ✓ Low retention



Sales unit :

768 Tips
