

TECHNICAL DATA SHEET

Part Number	041634
Description	Thermo Scientific™ ART™ Barrier Extended Length Pipette Tips



Benefits :

- ✓ Unlike standard filter tips, Thermo Scientific™ ART™ Barrier Extended Length Pipette Tips in Hinged Racks seal when exposed to potential contaminants. The pipette tips from the reload inserts fit securely into the environmentally friendly, reusable hinged racks and the covers may be used as hinged or lift-off style. The locking clasp prevents contamination by keeping remaining tips sterile and offers an easy, one-handed operation. ART hinged racks are available in low retention, micropoint, extended length, and other specialty styles.
- ✓ Stackable hinged racks save space in your lab
- ✓ Sterile tips are certified free from RNase, DNase and pyrogen
- ✓ Length : 9.09 cm
- ✓ Volume : 200 µL



Sales unit :

768 Tips