

Filter tips, 2 - 20 µl, PP, CE-IVD

Cleanroom quality for contamination-free analyses.

Perfect fit for the Transferpette® micropipette and many other pipettes.

The tip system for all requirements.

Pipette tips and filter tips from BRAND offer reliable quality for precise and demanding analyses. The raw materials are free from DiHEMDA und Oleamid. Tips are manufactured at BRAND under monitored cleanroom conditions and are automatically packaged. No parting agents or demolding aids are used. All racked tips and filter tips up to 1000 µl are free of DNA (< 40 fg), RNase (< 8.6 fg), endotoxins (< 1 pg) and ATP (< 1 fg).

BRAND offers pipette tips for volumes from 0.1 µl to 10 ml and filter tips for volumes from 0.1 µl to 1000 µl. They are tested for BRAND pipettes and many pipette from other manufacturers. See Selection Guide for more information.

BRAND pipette tips and filter tips are autoclavable at 121 °C (2 bar) acc. DIN EN 285 and CE-marked according to IVD-Directive 98/79 EC

Sterile tips are manufactured under certified BIO-CERT® quality. They are sterilized with beta radiation according to ISO 11 137 and the AAMI Guidelines.

Non-self-sealing filter tips from BRAND have a PE filter that is free from chemical additives. Permeability is controlled by the distinct combination of pore size and filter length, so that no aerosols can reach the pipette shaft.



| Cat.-No. | Sterility | Type of packaging | Volume range | Color-Code tip-tray |
|----------|-------------|-------------------|--------------|---------------------|
| 732608 | non-sterile | TipRack | 2 - 20 µl | yellow |
| 732628 | BIO-CERT® | TipRack | 2 - 20 µl | yellow |
| 732708 | non-sterile | TipBox | 2 - 20 µl | yellow |
| 732728 | BIO-CERT® | TipBox | 2 - 20 µl | yellow |