



## TECHNICAL DATA SHEET

### **PRODUCT DESCRIPTION**

<b>PRODUCT CODE</b>	162C
<b>GENERAL DESCRIPTION</b>	TAMPONE SECCO A/AL 1000PKG
<b>DETAILED DESCRIPTION</b>	10 boxes containing 100 tubes with aluminum stick with polyester tip and orange cap in bulk

### **PRODUCT COMPONENTS**

<b>APPLICATOR</b>	ALUMINIUM APPLICATOR	<b>DIMENSIONAL SPECIFICATIONS</b>
		STICK LENGTH: 150.0 mm $\pm$ 1.0 mm TOTAL LENGTH: 152.0 mm $\pm$ 1.0 mm STICK DIAMETER: 1.00 mm $\pm$ 0.10 mm TIP DIAMETER: 3.00 mm $\pm$ 0.50 mm TIP LENGTH: 15.00 mm $\pm$ 2.00 mm
<b>CAP</b>	CAP YELLOW	<b>DIMENSIONAL SPECIFICATIONS</b>
		CAP LENGTH: 32.15 mm $\pm$ 0.20 mm INSERTION DIAMETER: 9.30 mm $\pm$ 0.05 mm EXTERNAL DIAMETER: 12.2 mm $\pm$ 0.10 mm
<b>TRANSPORT TUBE</b>	TUBE	<b>DIMENSIONAL SPECIFICATIONS</b>
		TUBE LENGTH: 155.00 mm $\pm$ 1.00 mm INSERTION DIAMETER: 9.0 mm $\pm$ 0.5 mm EXTERNAL TUBE DIAMETER: 12.0 mm $\pm$ 0.5 mm

**TECHNICAL DATA SHEET****PACKAGING SPECIFICATION**

<b>INNER BOX QTY</b>	100
<b>OUTER BOX QTY</b>	1000
<b>PACKAGING</b>	BULK

**IMPORTANT INFO**

<b>STERILIZATION METHOD</b>	GAMMA IRRADIATION
<b>STORAGE TEMPERATURE</b>	5/25°C-DO NOT FREEZE
<b>SHELF LIFE</b>	36 MONTHS