



## TECHNICAL DATA SHEET

### **PRODUCT DESCRIPTION**

<b>PRODUCT CODE</b>	173CS01
<b>GENERAL DESCRIPTION</b>	SPREADERS A T B/SING
<b>DETAILED DESCRIPTION</b>	1 box containing 500 paper-film pouches containing one T-shape spreader

### **PRODUCT COMPONENTS**

<b>SPREADERS</b>	SPREDERS A T	<b>DIMENSIONAL SPECIFICATIONS</b>
		TOTAL LENGTH: 141.5 mm $\pm$ 0.5 mm DIAMETER: 4.0 mm $\pm$ 0.2 mm THICKNESS "T": 4.0 mm $\pm$ 0.2 mm LENGTH "T": 35.0 mm $\pm$ 0.5 mm

### **PACKAGING SPECIFICATION**

<b>PACKAGING</b>	PAPER-FILM POUCH
<b>BOX QTY</b>	500

### **IMPORTANT INFO**

<b>STERILIZATION METHOD</b>	GAMMA IRRADIATION
<b>STORAGE TEMPERATURE</b>	2/30°C – DO NOT FREEZE
<b>SHELF LIFE</b>	36 MONTHS