



## TECHNICAL DATA SHEET

### **PRODUCT DESCRIPTION**

<b>PRODUCT CODE</b>	178CS20
<b>GENERAL DESCRIPTION</b>	ANSE 1UL SOFT B/20PZ 4000PKG
<b>DETAILED DESCRIPTION</b>	4 boxes each containing 50 paper-film pouches each containing 20 1ul-loops

### **PRODUCT COMPONENTS**

		<b>DIMENSIONAL SPECIFICATIONS</b>
<b>LOOP</b>	ANSA 1UL SOFT LIGHT GREEN	TOTAL LENGTH: 199.00 mm $\pm$ 0,50 mm EXTERNAL DIAMETER RING: 3.00 mm $\pm$ 0.10 mm INTERNAL DIAMETER RING: 0.80 mm $\pm$ 0.10 mm THICKNESS RING: 1.10 mm $\pm$ 0.10 mm

### **PACKAGING SPECIFICATION**

<b>INNER BOX QTY</b>	1000
<b>OUTER BOX QTY</b>	4000
<b>PACKAGING</b>	PAPER-FILM POUCH
<b>INTERNAL ENVELOPE QTY</b>	20

### **IMPORTANT INFO**

<b>STERILIZATION METHOD</b>	GAMMA IRRADIATION
<b>STORAGE TEMPERATURE</b>	2/30°C – DO NOT FREEZE
<b>SHELF LIFE</b>	36 MONTHS