



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	179CS20
GENERAL DESCRIPTION	ANSE 10UL SOFT B/20PZ 4000PKG
DETAILED DESCRIPTION	4 boxes each containing 50 paper-film pouches each containing 20 light blue "10ul" loops

PRODUCT COMPONENTS

		DIMENSIONAL SPECIFICATIONS
LOOP	ANSA 10UL SOFT LIGHT BLUE	TOTAL LENGTH: 201.00 mm \pm 0.50 mm EXTERNAL DIAMETER RING: 6.50 mm \pm 0.10 mm INTERNAL DIAMETER RING: 3.90 mm \pm 0.10 mm THICKNESS RING: 1.20 mm \pm 0.10 mm

PACKAGING SPECIFICATION

INNER BOX QTY	1000
OUTER BOX QTY	4000
PACKAGING	PAPER-FILM POUCH
ENVELOPE QTY	20

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.01-13/10/2015) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.