



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	176CS20
GENERAL DESCRIPTION	NEEDLES B/20PCS 4000PKG
DETAILED DESCRIPTION	4 boxes each containing 50 paper-film pouches each containing 20 needles

PRODUCT COMPONENTS

SPREADERS	Needle	DIMENSIONAL SPECIFICATIONS
		TOTAL LENGTH: 199.00 mm \pm 0.50 mm TIP DIAMETER: 1.50 mm \pm 0.10 mm

PACKAGING SPECIFICATION

INNER BOX QTY	1000
OUTER BOX QTY	4000
PACKAGING	PAPER-FILM POUCH
ENVELOPE QTY	20

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS