



## TECHNICAL DATA SHEET

### PRODUCT DESCRIPTION

<b>PRODUCT CODE</b>	312C25T
<b>GENERAL DESCRIPTION</b>	LOOPS 4MM S/M T/25PZ MET.
<b>DETAILED DESCRIPTION</b>	1 tube with conical bottom and green screw cap containing 25 "4 mm" loops

### PRODUCT COMPONENTS

<b>LOOP</b>	LOOP 4 mm	<b>DIMENSIONAL SPECIFICATIONS</b>
		EXTERNAL LOOP DIAMETER: 4.05 mm $\pm$ 0.10 mm WIRE DIAMETER: 0.74 mm $\pm$ 0.05 mm TOTAL LENGTH: 75.0 mm $\pm$ 2.0 mm
<b>SAMPLING TUBE</b>	CONICAL TUBE	<b>DIMENSIONAL SPECIFICATIONS</b>
		TUBE LENGTH: 101.5 mm $\pm$ 0.2 mm TUBE DIAMETER: 16.30 mm $\pm$ 0.20 mm INSIDE DIAMETER: 14.1 mm $\pm$ 0.2 mm

This bulletin (Rev.02-07/11/2016) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



## TECHNICAL DATA SHEET

<b>CAP</b>	GREEN CAP	<b>DIMENSIONAL SPECIFICATIONS</b>
		CAP DIAMETER: 21.50 mm $\pm$ 0.20 mm CAP LENGTH: 16.30 mm $\pm$ 0.20 mm

### **PACKAGING SPECIFICATION**

<b>TUBE QUANTITY</b>	25
----------------------	----

### **IMPORTANT INFO**

<b>STERILIZATION METHOD</b>	NOT STERILE
<b>STORAGE TEMPERATURE</b>	2/30°C – DO NOT FREEZE
<b>SHELF LIFE</b>	60 MONTHS

This bulletin (Rev.02-07/11/2016) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.