



## TECHNICAL DATA SHEET

### PRODUCT DESCRIPTION

<b>PRODUCT CODE</b>	204C
<b>GENERAL DESCRIPTION</b>	PIP.PASTEUR 1,5ML 23CM 2400PKG
<b>DETAILED DESCRIPTION</b>	6 boxes containing 400 "1,5 ml" long Pasteur pipettes in bulk

### PRODUCT COMPONENTS

		<b>DIMENSIONAL SPECIFICATIONS</b>
<b>PIPET</b>	LONG PIPETTE 1,5 ML	DIAMETER BULB: 13.20 mm $\pm$ 0.30 mm TOTAL LENGTH: 227.00 $\pm$ 1.50 STEM DIAMETER: 4.40 mm $\pm$ 0.30 mm DROPS PER ML: 25 CAPACITY: 6 ml BULB DRAW: 3,2 ml GRADUATED UP TO: 1.5 ml GRAD.INTVL.: 0.5 ml

### PACKAGING SPECIFICATION

<b>INNER BOX QTY</b>	400
<b>OUTER BOX QTY</b>	2400
<b>PACKAGING</b>	BULK

### IMPORTANT INFO

<b>STERILIZATION METHOD</b>	NOT STERILE
<b>STORAGE TEMPERATURE</b>	2/30°C – DO NOT FREEZE
<b>SHELF LIFE</b>	60 MONTHS

This bulletin (Rev.01-02/09/2015) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.