



## TECHNICAL DATA SHEET

### **PRODUCT DESCRIPTION**

<b>PRODUCT CODE</b>	167KS01
<b>GENERAL DESCRIPTION</b>	PL.APPL.RAYON IND.WRAP.1000PCS
<b>DETAILED DESCRIPTION</b>	1 box containing 1000 paper-film pouches each containing 1 plastic applicator with rayon tip

### **PRODUCT COMPONENTS**

<b>APPLICATOR</b>	PLASTIC APPLICATOR	<b>DIMENSIONAL SPECIFICATIONS</b>
		STICK LENGTH: 144.0 mm $\pm$ 1.0 mm PLASTIC TIP DIAMETER: 2.50 mm $\pm$ 0.20 mm FIBRE TIP LENGTH: 15.0 mm $\pm$ 2.0 mm FIBRE TIP DIAMETER: 4.50 mm $\pm$ 0.50 mm TOTAL LENGTH: 147.0 mm $\pm$ 1.0 mm

### **PACKAGING SPECIFICATION**

<b>PACKAGING</b>	PAPER-FILM POUCH
<b>BOX QTY</b>	1000

### **IMPORTANT INFO**

<b>STERILIZATION METHOD</b>	GAMMA IRRADIATION
<b>STORAGE TEMPERATURE</b>	2/30°C – DO NOT FREEZE
<b>SHELF LIFE</b>	36 MONTHS