



## TECHNICAL DATA SHEET

### **PRODUCT DESCRIPTION**

<b>PRODUCT CODE</b>	948C
<b>GENERAL DESCRIPTION</b>	SRK 4ML UNS BULK 160MM TUBE
<b>DETAILED DESCRIPTION</b>	6 boxes containing 50 tubes in bulk with 4ml of medium UNS and 1 plastic applicator and green screw cap.

### **PRODUCT COMPONENTS**

<b>MEDIUM</b>	<b>COMPOSITION</b>
UNS MLS	Peptone Sodium Chloride Disodium phosphate Monopotassium phosphate Tween 80 Lecithin Sodium thiosulfate Sodium Thioglicollate Sodium Bicarbonate Potassium chloride

<b>APPLICATOR</b>	<b>PLASTIC APPLICATOR</b>	<b>DIMENSIONAL SPECIFICATIONS</b>
		STICK LENGTH: 144.0 mm $\pm$ 1.0 mm PLASTIC TIP DIAMETER: 2.50 mm $\pm$ 0.20 mm FIBRE TIP LENGTH: 15.0 mm $\pm$ 2.0 mm FIBRE TIP DIAMETER: 4.50 mm $\pm$ 0.50 mm TOTAL LENGTH: 147.0 mm $\pm$ 1.0 mm
<b>SAMPLING TUBE</b>	<b>TUBE</b>	<b>DIMENSIONAL SPECIFICATIONS</b>
		TOTAL LENGTH: 158.50 mm $\pm$ 0.20 mm TUBE DIAMETER: 16.00 mm $\pm$ 0.20 mm INTERNAL DIAMETER: 14.50 mm $\pm$ 0.20 mm

This bulletin (Rev.00-22/12/2014) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



## TECHNICAL DATA SHEET

CAP	GREEN CAP	DIMENSIONAL SPECIFICATIONS
		CAP DIAMETER: 21.50 mm $\pm$ 0.20 mm CAP LENGTH: 16.30 mm $\pm$ 0.20 mm

### **PACKAGING SPECIFICATION**

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	BULK

### **IMPORTANT INFO**

STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	12 MONTHS
MEDIUM pH	7.3 $\pm$ 0.5

This bulletin (Rev.00-22/12/2014) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.