

TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	940C
GENERAL DESCRIPTION	SRK SPONGE DEVICE 300PKG
DETAILED DESCRIPTION	6 boxes containing 50 paper-film pouches each containing a tube with a green screw cap and 10 ml of medium SRK with a blue spatula with sponge and blue cap.

PRODUCT COMPONENTS

MEDIUM		COMPOSITION
SRK		Sodium Chloride Potassium Chloride Calcium Chloride Lecithin Sodium thiosulphate Sodium Thyoglycollate Sodium Pyruvate Sodium Bicarbonate Polysorbate 80
SAMPLING TUBE	TUBE WITH ROUNDED BOTTOM	DIMENSIONAL SPECIFICATIONS
		TUBE LENGTH: 101.80 mm \pm 0.20 mm TUBE DIAMETER: 16.10 mm \pm 0.20 mm INTERNAL DIAMETER TUBE: 14.10 mm \pm 0.20 mm
CAP	GREEN CAP	DIMENSIONAL SPECIFICATIONS
		CAP DIAMETER: 21.50 mm \pm 0.20 mm CAP LENGTH: 16.30 mm \pm 0.20 mm



TECHNICAL DATA SHEET

SPONGE	SPONGE	DIMENSIONAL SPECIFICATIONS
		LENGTH OF SPONGE: 40.0 mm \pm 2.0 mm HEIGHT: 10 mm \pm 0.5 mm THICKNESS: 5.00 mm \pm 0.50 mm
CAP	CAP DUAL BLUE HOLE C\RECEPTACLE	DIMENSIONAL SPECIFICATIONS
		HEIGHT CAP: 16.60 mm \pm 0.20 mm CAP DIAMETER: 21.10 mm \pm 0.20 mm
APPLICATOR	BLUE SMALL SPATULA	DIMENSIONAL SPECIFICATIONS
		TOTAL LENGTH: 86.20 mm \pm 1.00 mm THICKNESS SPATULA: 0.80 mm \pm 0.10 mm WIDTH SPATULA: 6.90 mm \pm 0.20 mm

PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	PAPER-FILM POUCH

**TECHNICAL DATA SHEET****IMPORTANT INFO**

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	24 MONTHS
MEDIUM pH	7.0 ± 0.5