



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	920C
GENERAL DESCRIPTION	SRK 4ML LETHEEN 300PKG
DETAILED DESCRIPTION	6 boxes each containing 50 tubes with red screw cap containing a plastic applicator with rayon tip and 4 ml of medium LETHEEN BROTH

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
Lethen Broth	Peptone Sodium chloride Lecithin Tween 80 Sodio fosfato Beef extract

SAMPLING TUBE	TUBE WITH ROUNDED BOTTOM	DIMENSIONAL SPECIFICATIONS
		TUBE LENGTH: 101.80 mm \pm 0.20 mm TUBE DIAMETER: 16.10 mm \pm 0.20 mm INTERNAL DIAMETER TUBE: 14.10 mm \pm 0.20 mm
APPLICATOR	PLASTIC APPLICATOR	DIMENSIONAL SPECIFICATIONS
		STICK LENGTH: 83.30 mm \pm 0.50 mm DIAMETER HANDLE PART: 2.50 mm \pm 0.20 mm TIP LENGTH: 16.50 mm \pm 2.00 mm FIBRE TIP DIAMETER: 4.50 mm \pm 0.50 mm TOTAL LENGHT: 85.50 mm \pm 2.00 mm

This bulletin (Rev.03-02/04/2015) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

CAP	RED CAP	DIMENSIONAL SPECIFICATIONS
		CAP DIAMETER: 21.50 mm \pm 0.20 mm CAP LENGTH: 16.30 mm \pm 0.20 mm

PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	BULK

IMPORTANT INFO

STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	12 MONTHS
MEDIUM pH	7.0 \pm 0.5

This bulletin (Rev.03-02/04/2015) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.