



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	167CS01
GENERAL DESCRIPTION	DUAL SWABS PL.APP.IND.WRAPPED
DETAILED DESCRIPTION	10 boxes each containing 50 paper-film pouches each containing 1 tube with white cap and 2 plastic applicators

PRODUCT COMPONENTS

APPLICATOR	PLASTIC APPLICATOR	DIMENSIONAL SPECIFICATIONS
		TOTAL LENGTH: 145.50 mm \pm 1.00 mm STICK LENGTH: 141.00 mm \pm 1.00 mm PLASTIC TIP DIAMETER: 2.50 mm \pm 0.20 mm FIBRE TIP LENGTH: 17.00 mm \pm 2.00 mm FIBRE TIP DIAMETER: 4.50 mm \pm 0.50 mm
CAP	DUAL WHITE CAP	DIMENSIONAL SPECIFICATIONS
		EXTERNAL DIAMETER: 16.00 mm \pm 0.10 mm CAP LENGTH: 31.00 mm \pm 0.20 mm INSERTION DIAMETER: 12.00 mm \pm 0.10 mm
TRANSPORT TUBE	DUAL TUBE	DIMENSIONAL SPECIFICATIONS
		INTERNAL DIAMETER: 11.00 mm \pm 0.20 mm TOTAL LENGTH: 153.50 mm \pm 1.00 mm

This bulletin (Rev.00-13/02/2014) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.

**TECHNICAL DATA SHEET****PACKAGING SPECIFICATION**

INNER BOX QTY	50
OUTER BOX QTY	500
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	36 MONTHS