



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	3102CS01
GENERAL DESCRIPTION	B.TREFF 2ML+SPIN BASKET B/01
DETAILED DESCRIPTION	6 boxes containing 100 bags paper-film containing 1 PIPE 2 ml DNA FREE and 1 SPIN basket each ones.

PRODUCT COMPONENTS

SAMPLING TUBE	SPIN BASKET	SPECIFICATIONS
		HEIGHT: 23.42 ± 0.10 mm DIAMETER: 13.05 mm ± 0.05 mm MATERIALS: POLYPROPYLENE
SAMPLING TUBE	BRAND-TREFF 2 ML HDNA FREE	SPECIFICATIONS
		TUBE HEIGHT: 40.0 mm ± 0.2 mm TUBE EXTERNAL DIAMETER: 10.0 mm ± 0.2 mm CAP HEIGHT: 6.0 mm MATERIALS: POLYPROPYLENE

PACKAGING SPECIFICATION

INNER BOX QTY	100
OUTER BOX QTY	600
PACKAGING	PAPER-FILM POUCH

This bulletin (Rev.01-07/09/2018) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

IMPORTANT INFO

STERILIZATION METHOD	ETO TREATMENT
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS