

# Product Specifications

## General Purpose

Key#	Samco Catalog#	Product Description	Length	Stem Diameter	Overall Capacity	Bulb Draw	Drops per mL	Quantity Per Pack	Quantity Per Case	Sterile Pipets Individual	20's
1	202	General Purpose Standard	15.6 cm	7.9 mm	7.7 mL	3.2 mL	20	500	5,000	202-1S	202-20S
2	204	General Purpose Large Bulb	15.1 cm	6.1 mm	8.0 mL	4.6 mL	20	400	4,000	204-1S	204-20S
3	273	General Purpose 3"	8.7 cm	4.8 mm	1.7 mL	0.9 mL	25	500	4,000	273-1S	273-20S
4	335	5" B/B-Pet® (Blood Bank)	12.9 cm	6.6 mm	4.0 mL	1.9 mL	25	400	4,000	335-1S	335-20S
5	336	6" B/B Pet® (Blood Bank)	15.4 cm	6.6 mm	5.0 mL	1.9 mL	23	500	5,000	336-1S	336-20S
6	691	62 µL Large Aperture Pipet	13.7 cm	6.1 mm	9.3 mL	4.8 mL	—	250	2,500	691-1S	691-20S

## Graduated

Key#	Samco Catalog#	Product Description	Length	Stem Diameter	Overall Capacity	Bulb Draw	Drops per mL	Quantity Per Pack	Quantity Per Case	Sterile Pipets Individual	20's
7	212	Graduated 1 mL Small Bulb	14.0 cm	4.7 mm	3.9 mL	2.2 mL	25	500	5,000	212-1S	212-20S
8	222	Graduated 1 mL Large Bulb	15.5 cm	4.6 mm	5.8 mL	3.4 mL	25	500	5,000	222-1S	222-20S
9	225	Graduated 3 mL Large Bulb	15.5 cm	7.9 mm	7.5 mL	3.2 mL	20	500	5,000	225-1S	225-20S
10	274	Graduated 2 mL (Blood Bank)	15.2 cm	7.1 mm	4.8 mL	1.8 mL	24	500	5,000	274-1S	274-20S
11	282	Graduated 0.3 mL (Pediatric)	11.6 cm	3.0 mm	1.5 mL	0.8 mL	21	500	4,000	282-1S	282-20S

## Fine Tip

Key#	Samco Catalog#	Product Description	Length	Stem Diameter	Overall Capacity	Bulb Draw	Drops per mL	Quantity Per Pack	Quantity Per Case	Sterile Pipets Individual	20's
12	232	Fine Tip Standard	14.7 cm	4.6 mm	5.8 mL	3.3 mL	55	500	5,000	232-1S	232-20S
13	234	Fine Tip Large Bulb	14.7 cm	6.1 mm	8.7 mL	4.8 mL	45	400	4,000	234-1S	234-20S
14	231	Extended Fine Tip Small Bulb	10.4 cm	3.0 mm	1.5 mL	1.0 mL	50	400	3,200	231-1S	231-20S
15	233	Extended Fine Tip Large Bulb	15.3 cm	3.0 mm	5.0 mL	3.3 mL	50	250	2,000	233-1S	233-20S
16	235	Extended Fine Tip Standard	15.7 cm	4.6 mm	5.8 mL	3.5 mL	50	500	5,000	235-1S	235-20S
17	243	Fine Tip Small	5.1 cm	1.19 mm	1.3 mL	1.0 mL	45	500	4,000	—	

Custom design options available to meet your exact specifications