

TECHNICAL DATA SHEET

Part Number	043225
Description	Transfer pipettes



Benefits :

Single-use transfer pipettes are suitable for the vast majority manipulations of laboratory liquids. These pipettes are made of polyethylene and are therefore unbreakable. The bulb integrated into the pipette contributes to a greater safety of use. The transfer pipettes meet both the requirements of the work of routine than more specific work. Non-toxic polyethylene has a low affinity for proteins.

Technical specifications :

- ✓ Total volume : 7.5ml
- ✓ Bulb volume : 3.2 ml
- ✓ Length : 155 mm
- ✓ Drops per ml: 20

Sales unit :

500 Pipettes
