



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	224C
GENERAL DESCRIPTION	PIP.PAST FINE/STD BULB 3000PKG
DETAILED DESCRIPTION	1 box contain 3000 pipettes organized in 6 box each containing 500 pipettes.

PRODUCT COMPONENTS

PIPET	PIPETTE FINE TIP	SPECIFICATIONS
		TOTAL LENGTH: 146.0 mm \pm 1.0 mm STEM DIAMETER: 5.00 mm \pm 0.25 mm DIAMETER BULB: 13.50 mm \pm 0.30 mm PIPETTE MATERIALS: POLYETHYLENE

PACKAGING SPECIFICATION

INNER BOX QTY	500
OUTER BOX QTY	3000
PACKAGING	BULK

IMPORTANT INFO

SHELF LIFE	60 MONTHS
-------------------	-----------

This bulletin (Rev.00-26/09/2018) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.