

## TECHNICAL DATA SHEET

---

<b>Part Number</b>	043336S
<b>Description</b>	Samco™ Transfer Pipettes



---

### Benefits :

- ✓ Prevent cross-contamination with Thermo Scientific™ Samco™ General-Purpose Transfer Pipettes, which have a built-in bulb that eliminates the need for plugging, packaging and autoclaving.
- ✓ Made of shatterproof, nontoxic polyethylene.
- ✓ Their low-affinity surface reduces loss of cells and valuable proteins due to binding.
- ✓ They can be frozen in liquid nitrogen, or heat sealed for storage.
- ✓ Sterile
- ✓ Volume : 4.9 ml
- ✓ Bulb Draw : 1.8 ml
- ✓ Sterile

### Sales unit :

500 Pipettes

---