


Data Sheet 215006

This document declares that the 96 –well Krystal Black Assay Plate, Tissue Culture treated, Pack of 100 meets the following specification and drawing.

- Supplied lidded & sterile in individual bags
- SBS standard plate format
- Thickness of bottom 0.75mm
- Well volume of 350ul
- Constructed from ultra-pure grade Polystyrene
- Flat, optically clear base/bottom
- Assembled with ultra-sonic welding



Signed:

Katie Williams
Customer Service Co-ordinator
KWilliams@porvair-sciences.com

DOMINIQUE DUTSCHER SAS