

4titude[®] Sealing Solutions

4titude[®] Sealing Solutions

4titude[®] offers the widest range of plate sealing solutions available on the market. You can choose between sealing with strip caps, mats, lids, adhesive seals in strip or plate format, and heat seals in flexible formats up to plate size. The choice of an optimised sealing solution is particularly important for (q)PCR because thermal cycling can be associated with evaporation of reaction reagents.

Please see below for a brief overview of our sealing solutions for PCR plates and strips and refer to our webpage www.4ti.co.uk/seal or the **4titude[®] Sealing Solutions Product Overview brochure** for more details.

Strips of Sealing Caps

- Available in strips of 8 or 12, either domed or flat optical caps
- Universal fit for all 4titude[®] PCR tube strips and PCR plates
- CrystalStrips[™] made of a special polymer with improved optical properties leading to high transmission rates, ideally suited for small samples with low signal intensity

Code	Details	Strips/Case
4ti-0751	Strips of 8 flat optical caps	300
4ti-0783	as above	125
4ti-0752	Strips of 8 domed caps	300
4ti-0782	as above	125
4ti-0755	CrystalStrip [™] , crystal clear strips of 8 flat optical caps	300
4ti-0755/125	as above	125
4ti-0788	Strips of 12 flat optical caps	200