


## FrameStar® 384 Well Skirted PCR Plate



### Features


- Compatible with the majority of 384 well block PCR, qPCR and sequencing instruments
- Available with frosted wells for increased qPCR signal intensities

### FrameStar® 384 Well Skirted PCR Plate

PP wells, PC frame, cut corner A24,  
working volume: <30 µl, total well capacity: 55 µl

Code	Details			Plates/Case	
4ti-0384		purple frame		clear wells	50
4ti-0384/B		blue frame		clear wells	50
4ti-0384/C		clear frame		clear wells	50
4ti-0384/G		green frame		clear wells	50
4ti-0384/R		red frame		clear wells	50
4ti-0384/X		black frame		clear wells	50
4ti-0385		black frame		white wells	50
4ti-0386		black frame		black wells	50
4ti-0387		clear frame		frosted wells	50

## FrameStar® 384 Well Skirted PCR Plate, Roche Style



### Features

- Designed for use on the Roche LightCycler® 480 with 384 well block
- Available as a combi pack with qPCR adhesive seals

### FrameStar® 384 Well Skirted PCR Plate, Roche Style

PP wells, PC frame, cut corner A24/P24,  
working volume: <30 µl, total well capacity: 55 µl

Code	Details			Plates/Case	
4ti-0380/C		clear frame		clear wells	50
4ti-0381		clear frame		white wells	50

Code Combi Pack		Plates and Seals/Case
4ti-0382	4ti-0381 + qPCR Seal (4ti-0560)	50
4ti-0383	4ti-0380/C + qPCR Seal (4ti-0560)	50