

[🔗](#) > KIMBLE® KONTES® DUALL® Tissue Grinder, with glass pestle and tube, Size 21, 3 mL



KIMBLE® KONTES® DUALL® Tissue Grinder, with glass pestle and tube, Size 21, 3 mL

- Manufactured from borosilicate glass 3.3
- This ground glass model is ideal for connective tissue such as muscle, heart and lung
- Construction is strong enough to allow the pestle to be motor-driven
- Pestles are designed to be used with 788000 stirrer adapters
- Replacement components are available and completely interchangeable
- Replacement tube: 885452 serie; replacement pestle: 885451 serie
- Clearance between pestles and tubes is 0.004 " to 0.006 "