

Specification Sheet



REUSABLE SCREW THREAD CULTURE TUBES

KIMAX® culture tube made from tubing with uniform wall thickness for maximum heat transfer and chemical resistance.

- Provided with phenolic screw-thread caps, unattached
- Series 45066 is designed from ASTM E982, Type VI, Class A requirements
- Manufactured from 33 expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements

Part Number	45066-16100
Name	Reusable Screw Thread Culture Tubes with Unattached Caps
Capacity to Neck (mL)	12
OD (mm)	16
Length (mm)	100
GPI Finish	15-415
Liner Material	Rubber
Shelf-Pack Qty	72
Case Qty	288