

Specification Sheet



REUSABLE SCREW THREAD CULTURE TUBES

KIMAX® culture tube made from tubing with uniform wall thickness for maximum heat transfer and chemical resistance.

- The 20 x 125 mm size of these tubes is suitable for the determination of insulin according to the method of Alving, Rubin and Millor, Journal of Biological Chemistry, 127, 3 (March 1939)
- Provided with phenolic screw-thread caps, unattached
- Series 45066 is designed from ASTM E982, Type VI, Class A requirements
- Manufactured from 33 expansion, low extractable borosilicate glass conforming to USP Type I and ASTM E438, Type I, Class A requirements

Part Number	45066-20125
Name	Reusable Screw Thread Culture Tubes with Unattached Caps
Capacity to Neck (mL)	25
OD (mm)	20
Length (mm)	125
GPI Finish	18-415
Liner Material	Rubber
Shelf-Pack Qty	48
Case Qty	192