



For Information only

ITEM	QTY	TITLE
2	1	CAP: BLACK PHENOLIC   LINER: .050" RUBBER
1	1	VIAL: TYPE 1 CLEAR BOROSILICATE GLASS, 33 EXP. PER ASTM E438
MATERIAL: DESCRIPTION AND SPECIFICATIONS		

12X35 CLEAR VIAL W/PHEN CAP



ITEM NUMBER	60811D-12	PARAMOUNT NUMBER	
OLD MATERIAL NUMBER			

