

Tubes shown in actual size.



Pirouet™ 0.5, 1.5, and 2 mL Microcentrifuge Tubes

- Two cap styles to suit customer preference:
 - Boil-proof caps are ideal for rigorous high-temp applications and tight storage conditions.
 - Large, boil-resistant caps are easy to open & close, with a large labeling surface.
- Crystal-clear polymer facilitates visualization of samples.
- Rated to 20,000 x G.
- Graduations at every 100 µL for easy volume estimation.
- Available with proprietary NoStick® resin to ensure maximum sample recovery.
- Temperature range: -80°C to +121°C.
- Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free. Sterile product is certified pyrogen-free.

Graduated Microcentrifuge Tubes

Item #	Max Fill	Description	Colors	Packaging
1110-00	0.65 mL	0.5 mL Boil-Proof Cap	-00, -01, -02, -03, -04, -05, -06, -07, -09, -29	1000/unit, 10 units/case
1210-00	1.7 mL	1.5 mL Boil-Proof Cap	-00, -01, -02, -03, -04, -05, -06, -09	500/unit, 10 units/case
1260-00	1.7 mL	1.5 mL Boil-Resistant, Large Cap	-00, -01, -02, -03, -04, -05, -06, -09	500/unit, 10 units/case
1310-00	2.0 mL	2.0 mL Boil-Proof Cap	-00, -01, -02, -03, -04, -05, -06, -07, -09	500/unit, 10 units/case

Graduated Microcentrifuge Tubes, Sterile

1110-50	0.65 mL	0.5 mL Boil-Proof Cap, Sterile	-00 Clear	500/unit, 10 units/case
1210-50	1.7 mL	1.5 mL Boil-Proof Cap, Sterile		250/unit, 10 units/case
1310-50	2.0 mL	2.0 mL Boil-Proof Cap, Sterile		250/unit, 10 units/case

NoStick®, Graduated Microcentrifuge Tubes



1110-10	0.65 mL	0.5 mL NoStick, Boil-Proof Cap	-00 Clear	500/unit, 10 units/case
1210-10	1.7 mL	1.5 mL NoStick, Boil-Proof Cap	-00 Clear	250/unit, 10 units/case
1310-10	2.0 mL	2.0 mL NoStick, Boil-Proof Cap	-00 Clear	250/unit, 10 units/case

Graduated Microcentrifuge Tubes without Caps

1270-00	1.7 mL	1.5 mL Capless	-00 Clear	500/unit, 10 units/case
1370-00	2.0 mL	2.0 mL Capless	-00 Clear	500/unit, 10 units/case

