



1210-00

1.5 mL Graduated Microcentrifuge Tubes with Boil-Proof Cap

DESCRIPTION

Product Description

Volume: 1.5 mL

Max Fill: 1.7 mL

Tube Description: Graduated Microcentrifuge Tubes

Cap Description: Boil-Proof Cap

Packaging: 500/unit, 10 units/case

1210-00S, 1210-00

- Boil-proof caps are ideal for rigorous high-temp applications and tight storage conditions.
- Crystal-clear polymer facilitates visualization of samples.
- Rated to 20,000 x G.
- Graduations at every 100 μ L for easy volume estimation.
- Available with proprietary NoStick®resin to ensure maximum sample recovery.
- Temperature range: -80°C to +121°C.
- For use in pharmaceutical industry, food industry, and molecular biology.
- Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free.

DOMINIQUE DUTSCHER SAS