


Flipper Racks have been a mainstay of laboratories around the world for more than 25 years. The innovative design of Flipper Racks provides researchers with greater flexibility on the bench without managing multiple racks. From the basic 4-Way Flipper rack design to the Sliplock and Revolver, researchers are sure to find the ideal Flipper Rack for their lab.

4-Way Flipper® Racks

- Small footprint is perfect for the benchtop or under the hood.
- Suitable for use in ultra-low temperature freezers (-90°).
- Available in standard or interlocking versions.
- Accommodates all popular tube sizes, from 0.5 mL to 50 mL, and 20-25 mm.
- Rugged polypropylene construction.



Press the interlocking racks together to form a larger, custom rack.

5410-29

Item #	Description	Tube Options	Colors	Packaging
5410-29	4-Way Interlock Flipper 171 x 91 x 54 mm	Side 1: 32 x 1.5-2.0 mL Side 2: 32 x 0.5 mL Side 3: 12 x 5 or 15 mL Side 4: 4 x 50 mL	Fluorescent: Blue, Green, Yellow, Orange, Pink	5/unit, 4 units/case
5415-59	4-Way Flipper 5 mL 171 x 91 x 54 mm	Side 1: 18 x 5 mL and 32 x 0.5 mL Side 2: 32 x 1.5-2.0 mL Side 3: 12 x 5 or 15 mL Side 4: 4 x 50 mL	Standard: Blue, Green, Red, Yellow	4/unit, 4 units/case
5420-29	4-Way Flipper 50 mL 8 x 50 mL Tubes 198 x 114 x 95 mm	Side 1: 8 x 50 mL Side 2: 20 x 20 mm Side 3: 24 x 5 or 15 mL Side 4: 40 x 1.5-2.0 mL	Fluorescent: Blue, Green, Yellow, Orange, Pink	5/unit, 2 units/case
5430-39	4-Way Flipper 25 mm 12 x 25 mm Tubes 198 x 114 x 95 mm	Side 1: 12 x 25 mm Side 2: 20 x 20 mm Side 3: 40 x 1.5-2.0 mL Side 4: 24 x 5 or 15 mL	Standard: Blue, Green, Red, Orange, Violet	5 /unit, 2 units/case



5410-29



5420-29

4-Way Flipper Racks with Tray

- Multi-tube compatibility.
- Four 4-way racks in a convenient tray.
- Each Flipper can be independently oriented.
- 4-way rack can be used either in the tray or independently.
- Standard or fluorescent colors.



5400-39

Item #	Dimensions	Tube Options	Colors	Packaging
5400-39	4-Way Flipper with Tray 243 x 102 x 41 mm	12 x 0.5 mL 9 x 1.5-2.0 mL 6 x 15 mL 2 x 50 mL	Standard: Green, Blue, Red, Yellow	1 tray/unit, 5 units/case
5400-29	243 x 102 x 41 mm	12 x 0.5 mL 9 x 1.5-2.0 mL 6 x 15 mL 2 x 50 mL	Fluorescent: Green, Yellow, Orange, Pink	1 tray/unit, 5 units/case



5400-29



Made in USA