




EUROTUBO® Pipetting bulb

This one-handed design is the simplest safety pipette filler to use available. Manufactured in natural orange rubber. Approx. drawn capacity: 25 ml. Single hand use, only two operating points. Evacuate via the automatic valve. Standard model, accommodates all pipettes. Ability to clean inside of bulb by removing patented valve and rinsing out.



1. Evacuate the air by pressing the bulb, as the drawing indicates
2. Intake by pressing on point B (↑)
3. Drain the liquid by pressing on point A (↓)
4. Blow out the pipette by pressing as the drawing indicates (point C).

code	description	bag quantity	bag weight	bag volume
19200	pipetting bulb	1	0.05	0.0004



Pipetting bulb

Made of rubber.
Used to avoid mouth pipetting and contamination risk.
Can be opened, cleaned and **autoclaved**.
Ideal for Wintrobe and Westergren pipettes.



code	description	bag quantity	bag weight	bag volume
19201	red pipetting bulb	1	0.04	0.0002



Pipette pumps

Several models for various pipette volumes.
Designed for fast and efficient pipetting with simple, one handed operation.
Pipettes fit smoothly into collar.
Rotate the knurled thumb wheel on the side for precision filling or dispensing, and press the fast release lever for quick emptying.
Easy to use and easily disassembled for cleaning.
Sizes are colour coded.
Pipette pumps resist acids and alkalis.



code	description	colour	bag quantity	bag weight	bag volume
W-100	up to 2 ml	blue	1	0.06	0.0002
W-110	up to 10 ml	green	1	0.06	0.0002
W-120	up to 25 ml	red	1	0.06	0.0002