



Bel-Art Products

Wayne, New Jersey 07470

Tel: 973-694-0500 * Fax: 973-694-7199

Internet Address: <http://www.bel-art.com>

Technical Data Sheet

Product #: 13595-0000

Description: Electrophoresis Gel Plate Drying Rack

Material: Aluminum, Blue Epoxy Coating, Neoprene Rubber, PVC

Dimensions: 7 3/4" x 7" x 7"H

Other: Minimum Width of plates 6 1/2"
7 slots are .31" (7.9mm) wide 3 slots are .465" (11.8mm) Distance
between slots .375" (9.5mm)