



661 Route 23 South, Wayne, New Jersey 07470
Tel: 973-694-0500 * Fax: 973-694-7199
Internet Address: <http://www.bel-art.com>

This is to verify that all items noted below are in conformance with the Purchase Order Contract, drawings, specifications, and other applicable documentation; that all process certifications, material certifications, are on file at this or another Bel-Art facility and are available for review by the customer and or Government Representatives upon request.

Product #: 18872-0000 **Mfg. Date:** NA **Lot #:** General

Description: Radioimmunoassay Tube Rack

Material: Polypropylene, 20% Talc White, FDA Grade

Dimensions: 355mm x 260mm x 89mm, (120 places = 12 x 10)
Average hole diameter 21.6mm (holds tubes 16-20mm)

Other: Steam Autoclavable @121°C for 20 minutes 15 p.s.i.g.,
Resin meets FDA regulation 21 CFR 177.1520. Based on the process and materials involved in the manufacturing of this product, this product is considered to be BSE/TSE protein free. Based on the information provided to us this product does not contain or come in contact with lead, latex, or other contaminants. To the best of our knowledge based on information provided to us this Product does not contain or come in contact with any animal origin.

Manufactured: Bel-Art Products, Pocomoke City, MD USA

Reviewed and Approved By:

Signature: 

Printed Name: Carl Sisco

Title: Technical Support Manager

Date: 8/2/2018

SUBSIDIARIES: MADDAK, INC * APPLIED COATING, INC
ADDITIONAL PLANTS IN POCOMOKE, MD * KANSAS CITY, MISSOURI