



Wayne, New Jersey 07470  
Tel: 973-694-0500 \* Fax: 973-694-7199  
Internet Address: <http://www.bel-art.com>

## *Technical Data Sheet*

- Product #:** 13166-000
- Description:** Red Biohazard Bag  
(with biohazard warning label)
- Material:** High Molecular Weight, High Density Polyethylene
- Approximate Dimensions:** 22cm x 28cm  
8 1/2" x 11"  
.03mm (.0018") .72mil Thick. Max. Volume .43 gallons
- Other:** Printed with sterilization indicator patch  
Autoclavable @ 121° C (250° F), Recommended operating time 40  
minutes at sustained heat level (Cold Start)  
Melting Point for HMHD Bags 127°C (260° F)  
Dart-Impact-165 grams per ASTM D1709-01  
Tensile Strength ASTM D1922-00a Not Required  
Holder not included, Instructions in English, French, German,  
Spanish., Store in a cool place out of direct sunlight.  
Made in USA

