

FICHE TECHNIQUE

# Finntip 1000µl

## Référence

Référence	Conditionnement
9401030	1000 / sachet
9401010	10 000 / vrac
9401070	200 / boîte
9401110	10 x 96 / rack
9401113	10 x 96 / rack, stérile
9401115	10 x 96 / recharge
94052410	10 x 96 / recharge, filtre stérile

## Description du produit

- Finntip 1000 / Finntip Filtre 1000 µl, stérile

## CE

- Oui

## Utilisation

- Les pointes Finntip et les pipettes forment un système de transfert de liquide.

## Composition

- Matériel : Polypropylène (PP), filtre PE

## Caractéristiques

- Volume : 100-1000 µl
- Couleur : Transparent
- Graduation : Non
- Filtre : Non (sauf 94052410)
- Stérilité : Non (sauf 9401113 et 94052410, Certificat de pureté inclus dans le packaging)
- Longueur : 71 mm
- Largeur: diamètre orifice pointe extérieur : 1,3 mm, diamètre orifice au niveau de l'embout : 10 mm
- Poids net : 0,91 g, avec filtre 0,99 g

## Compatibilité pipette

- E.g. Finnpiquette™ F1/F2/ Novus 100-1000 µl

Visitez notre site pour plus d'information :

[thermofisher.com/finntip](http://thermofisher.com/finntip)

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific - Thermo Electron LED SAS - 10, rue duquay trouin - 44807 Saint-Herblain - Tél. 02 28 03 20 00

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.