

TECHNICAL DATA SHEET

Part Number	047911
Description	Tips 10 ml



Benefits :

- ✓ Finntip Flex 10 ml ext – extra long
- ✓ Non filtered
- ✓ Colour : natural
- ✓ Packaging : bulk

Certificate :

- ✓ Non sterile
- ✓ Autoclavable
- ✓ Certification : IS 9001, ISO 14001, ISO 13485
- ✓ Certified endotoxin, DNA, DNase and RNase free

Technical specifications :

- ✓ Volume : 10 ml
- ✓ Length : 26,7 cm

Sales unit :

100 Tips
