



## GripSeal® Caps With Molded O-Ring

**SKU:** 2005-00

**Colors:** -00 Clear, -01 Blue, -02 Green, -03 Orange, -04 Red, -05 Violet, -06 Yellow, -40 White

**Assorted Colors:** -09 includes Blue, Green, Orange, Red, Violet, Yellow

**Packaging:** 500/unit, 10 units/case

2005-00, 2005-00S, 2005-01, 2005-01S, 2005-02, 2005-02S, 2005-03, 2005-03S, 2005-04, 2005-04S, 2005-05, 2005-05S, 2005-06, 2005-06S, 2005-40, 2005-40S, 2005-09, 2005-09S

- Made from ultra-pure polypropylene (meets FDA 21 CFR 177.1520 & USP Class VI) with an inert TPE O-ring molded directly into the cap.
- Rubber ribs on cap allow for effortless opening and closing—even when frozen.
- Contrasting color stripes enable easy branding or organization.
- Easy  $\frac{3}{4}$  turn to open and close.
- Easy labeling surface.

**GripSeal®**

