

## WHEATON® Natural Leak Resistant Bottles

### Quality and precision leak resistant bottles for IVD kitting and storage

WHEATON leak resistant bottles provide a quality and precision design, with a superior double seal closure for your IVD primary packaging requirements. Options available include narrow or wide mouth, and amber or natural, to suit your application



- High or Low Density Polyethylene, or Polypropylene ■ Amber or natural color ■ Available in narrow or wide mouth ■ Superior double seal closure for leak resistance ■ No-drip pour lip for safe and clean pouring ■ Polypropylene cap attached

### Product Specification

Item Code	Bottle Colour	Bottle Material	Capacity ml	Cap size	Description	Diameter x height mm	Inner Pack Qty	Neck	Pack quantity
209665	Natural	PP	30	28-410	30 mL Star Btl, WM PP Nat, 28-410 Cap	33 x 69	12	Wide Mouth	72