

PRODUCT SPECIFICATION

Product Code: 224982

Product Description: 4 mL Sample Vial In Lab File, Amber , Rb Ln

Brand: Wheaton

Type of Glass: Type I Borosilicate

Product Specification:

Cap/ Closure Size	13-425
Capacity (ml)	4ml
Capacity (ml)	4ml
Diameter (mm)	15
Height (mm)	48

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of SciLabware Limited by:

Signed: M.Ashley
Quality Technician

Date: 12 Feb 2020