
WHEATON® LAB FILE® Sample Vials, 8mL Standard Vials With Caps Attached, Amber, 14B Rubber Septa



WHEATON® LAB FILE® Sample Vials, 8mL Standard Vials With Caps Attached, Amber, 14B Rubber Septa

- Clear vials are converted from Type I borosilicate glass tubing
- Amber vials are converted from amber borosilicate glass tubing that conforms to USP Type I requirements for light transmission to protect light sensitive products
- Sample vials with caps attached to vial help maintain cleanliness
- Lab File® with partitioned trays provides an easy way to inventory samples or to store empty vials

- Shorty Vials® feature a low profile for greater stability and less dead volume
- 14B rubber lined black phenolic screw cap for aqueous samples
- PTFE / 14B rubber lined black phenolic screw cap for organic samples
- Replacement caps can be purchased separately

Catalog Number: **224984**

DOMINIQUE DUTSCHER SAS