

TECHNICAL DATA SHEET

Part Number	048407
Description	WHEATON® Sample Vials, 40mL With Caps Attached



Benefits :

- ✓ Amber vials are converted from amber borosilicate glass tubing that conforms to USP Type I requirements for light transmission to protect light sensitive products
- ✓ Sample vials with caps attached to vial help maintain cleanliness
- ✓ Lab File® with partitioned trays provides an easy way to inventory samples or to store empty vials
- ✓ Shorty Vials® feature a low profile for greater stability and less dead volume
- ✓ 14B rubber lined black phenolic screw cap for aqueous samples
- ✓ PTFE / 14B rubber lined black phenolic screw cap for organic samples

Technical specifications :

- ✓ Dimensions (Diameter x height) : 28 x 98 mm
- ✓ Color : amber
- ✓ Capacity : 40 mL
- ✓ Thread : 24-400
- ✓ Material : Glass type: Amber Type 1A Borosilicate
- ✓ Cap Material : Black Phenolic
- ✓ Liner Material : PTFE / 14B Rubbe

Sales unit :

72 Sample vials