

TECHNICAL DATA SHEET

Part number	048445
Description	Glass scintillation vial



Description :

- ✓ Material : ASTM Type I Borosilicate
- ✓ Color : clear
- ✓ Finish Type Screw Thread
- ✓ Cap Material : Urea
- ✓ Thread Finish : 22-400
- ✓ Cap Liner Material : Metal Foil
- ✓ Cap/ Closure Size 22-400
- ✓ Capacity (mL or dr) : 20mL or 5.4dr
- ✓ Diameter : 28 mm
- ✓ Height : 61 mm

Piece by sales :

500 Vials
