

## Product data – LYPHAN boxes

**Product Name:** LYPHAN boxes L+ No. (e.g. L674; L666-8..)

**Country of origin:** Federal Republic of Germany

**Customs tariff No.:** 38220000

**SKU information:** 1 box = 200 strips

**Product description:**      **Product itself**                      **Product with plastic box**

Dimensions:                      85 mm x 11 mm                      80 mm x 96 mm x 18 mm

Paper thickness:                      0,27 - 0,30 mm

Weight of the strip/box:      0,14 g                                      ~ 65 g with 200 strips

Components:                      filter paper; indicator; cotton dyes

No minimum order:              1 box with 200 strips

No UPC or EAN; no UN number; no UNSPSC; no barcode

no dangerous product; SDS attached

**Storage temperature:**      room temperature

**Shelf life:**                      at least 12 month after receipt by consignee usable, if the strips stored under dry and dark conditions.

**Product picture:**

